# CALASEPT® EDTA

EDTA Solution 17%

CE  
0537

Magnum Dental OÜ,  
Aardla 13, 50112 Tartu, Estonia

**DIRECTA**

Distributed by

**DIRECTA AB**

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## 17% EDTA Root Canal Irrigant

### Indications

- Chemical and mechanical root canal preparation
- Removal of smear layer and inorganic substances
- Identification of root canals and canal openings
- Enlarging of root canals

### Composition:

17% EDTA in purified water

### Instructions for use

1. Apply Calasept EDTA at entrance of canal, file and ream canal in conventional manner.
2. Rinse alternately with sodium hypochlorite in order to remove pulp rests, necrotic tissue or dentine residues. Reapply Calasept EDTA as necessary until canal is enlarged to the desired size.
3. Rinse canal thoroughly so that all Calasept EDTA is removed from the canal.

### Contra indications:

Allergic reactions to edetic acid/sodium edetate.  
Wide open foramen apicale.

### Side effects:

No side-effects are known if properly applied. Irritation of periapical tissue is likely to occur if the material is spread over the apex.

### Interactions with other substances:

Unknown.

### Additional notes:

Keep out of the reach of children.  
For dental use only.

### Shelf life:

Two years under proper storage. This product must not be used after its shelf life has expired.

### Storage:

< 25°C. Store at room temperature. Always close the lid immediately after use.

**P305** IF IN EYES; **P351** Rinse cautiously with water for several minutes.  
**P338** Remove contact lenses, if present and easy to do. Continue rinsing. Refer to special instructions/Safety data sheet.

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